


## Our Scalers

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 18 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. of course, any manufacturer can make these claims, but we guarantee and certify ours.



- High-polish mirror finish tips.
- Manufactured within the guidelines of the FDA quality assurance system for good manufacturing practices.
- Highest grade German stainless steel, razor sharp, durable tips. Rockwell hardened to 60 RHc.
- 27 of the most popular designs from gracesys to McCalls to columbias, and many european patterns.
- All instruments are guaranteed against defects in material and workmanship.

- All instruments bear the ce mark, emblematic of the european community's highest standards.
- Customer satisfaction guarantee, too!
- Less than 20 grams ... Lightweight for less hand fatigue.
- Manufactured to the most exacting standards so the 50th instrument you buy is exactly like your first.
- You decide ... Why pay more for the same high-quality instruments?

## VALUED BUYERS

Compare our scalers and curettes with the best in the industry and you will immediately recognize we are equal in every way but one... price. Note the deep precision knurling for better tactile feel. The same high-grade stainless steel tips, the same light weight (under 20 grams) for less hand fatigue, the same professional craftsmanship, and surpassing the highest industry standards. Of course, any manufacturer can make these claims, but we guarantee and certify ours. We regularly manufacturing below articles in showing handles & pattern A to L for better grip & perfection. We are also fully equipped & capable to develop these instruments as per your own specific requirements & demands in materials, sizes, shapes and patterns.



A-Hollow Handle Ø 10 mm



B-Hollow Handle Ø 10 mm



C-Hollow Handle Ø 10 mm



D-Hollow Handle Ø 10 mm



E-Hollow Handle Ø 10 mm



F-Hollow Handle Ø 10 mm



G-Hollow Handle Ø 10 mm



H-Hollow Handle Ø 6.4 mm



I-Hollow Handle Ø 8 mm



J-Hollow Handle Ø 10 mm



K-Soild Handle Ø 6 mm



L-Soild Handle Ø 6 mm



- зеркал
- Espejos de la boca
- Mundspiegel
- дзеркал
- Miroirs à la bouche
- مرآيا الفم

## Plain



**DE-600** Fig. 3 Ø20mm  
**DE-601** Fig. 4 Ø22mm  
**DE-602** Fig. 5 Ø24mm  
**DE-603** Fig. 6 Ø26mm



## Cone-Socket Plain



**DE-604** Fig. 3 Ø20mm  
**DE-605** Fig. 4 Ø22mm  
**DE-606** Fig. 5 Ø24mm

## Magnifying



**DE-607** Fig. 3 Ø20 mm  
**DE-608** Fig. 4 Ø22 mm  
**DE-609** Fig. 5 Ø24 mm



## Cone-Socket Magnifying



**DE-610** Fig. 3 Ø20mm  
**DE-611** Fig. 4 Ø22mm  
**DE-612** Fig. 5 Ø24mm

## Rhodium



**DE-613** Fig. 3 Ø20mm  
**DE-614** Fig. 4 Ø22mm  
**DE-615** Fig. 5 Ø24 mm



## Cone-Socket Rhodium



**DE-616** Fig. 3 Ø20mm  
**DE-617** Fig. 4 Ø22mm  
**DE-618** Fig. 5 Ø24mm

## Front Surface



**DE-619** Fig. 3 Ø20mm  
**DE-620** Fig. 4 Ø22mm  
**DE-621** Fig. 5 Ø24mm  
**DE-622** Fig. 6 Ø26mm



## Cone-Socket Front Surface



**DE-623** Fig. 3 Ø20mm  
**DE-624** Fig. 4 Ø22mm  
**DE-625** Fig. 5 Ø24mm

## HANDLES FOR MOUTH MIRRORS

- Ручки для зеркал
- Manches pour miroirs à bouche
- Mundspiegelgriffe
- Ручки для дзеркал
- Mangos para espejos de boca
- مقابض لمرآيا الفم



**DE-626**  
Hollow Mirror  
Handle Round  
Ø 10 mm



**DE-627**  
Hollow Mirror  
Hexagonal  
Handle  
Ø 8 mm



**DE-628**  
Hollow  
Ø 8 mm



**DE-629**  
Hollow Mirror  
Hexagonal  
Handle  
Ø 10 mm



**DE-630**  
Hollow Mirror  
Handle  
Ø 8 mm



**DE-631**  
Hollow Mirror  
Handle  
Ø 10 mm



Cone-Socket Mirror Handle



**DE-632**  
Mirror Handle  
Solid  
Ø 6 mm



**DE-633**  
Mirror Handle  
Solid  
Ø 6 mm



**DE-634**  
Mirror Handle  
Solid Octagonal  
Ø 6 mm



**DE-635**  
Mirror Handle  
Stem  
with Rule  
35mm



**DE-636**  
Mirror Handle  
Solid  
Stem with  
Rule 35mm



**DE-637**  
Cone-Socket  
Solid Octagonal  
Ø 6 mm



**DE-638**  
Cone-socket  
Solid  
Ø 6 mm



**DE-639**  
Cone-socket  
Hexagonal  
Handle  
Ø 6 mm

- Зонды
- Sondes
- Zahnsonden

- Зонди
- Sondas
- المستكشفين

Examination kit Set of 3 Pieces



DE-640  
Probe

DE-641  
Mirror

DE-642  
College Forcep  
15cm



Examination Set of 4 Pieces



DE-643  
Exp,23/17-A

DE-644  
Exavators

DE-645  
Tweezers

DE-646  
Mouth Mirror

DE-647

First Aid Examination Set of 3 Pieces



DE-648

First Aid Examination Set of 3 Pieces



DE-649

First Aid Examination Set of 3 Pieces



DE-650

First Aid Examination Set of 3 Pieces



## TWEEZERS

- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten

- Пинцети
- Pinzas para algodón
- ملاقط



**DE-651**  
Surgical  
11.5cm



**DE-652**  
London College  
15cm



**DE-653**  
Semken Taylor  
Str.



**DE-654**  
Semken Taylor  
Cvd.



**DE-655**  
Meriam  
16cm



**DE-656**  
Dressing  
14cm



**DE-657**  
Semken Taylor  
1x2 Straight



**DE-658**  
Semken Taylor  
1x2 Curved

MEDSAL  
INTERNATIONAL



**DE-659**  
Gilles  
15cm 1x2



**DE-660**  
14cm 1x2



**DE-661**  
Perry  
13cm



**DE-662**  
Atraumatic  
15cm Cvd



**DE-663**  
Atraumatic  
15cm Str.



**DE-664**  
College with Lock  
15cm Cvd.



**DE-665**  
Gerald  
17.5cm



**DE-666**  
with lock  
15cm



**DE-667**  
Atraumatic  
17.5cm



# TWEEZERS

- Пинцеты
- Précelles - Brucelles
- Zahnpinzetten

- Пинцети
- Pinzas para algodón
- ملاقط



**DE-668**  
Adson-TC  
12cm



**DE-669**  
Adson-TC  
12cm 1x2



**DE-670**  
Adson-TC  
15cm 1x2



**DE-671**  
Adson  
12cm



**DE-672**  
Adson  
12cm 1x2



**DE-673**  
Adson  
15cm 1x2



**DE-674**  
Adson Brown  
12cm



**DE-675**  
Crane Kaplan  
15cm



**DE-676**  
Crane Kaplan  
15cm



**DE-677**  
Crane Kaplan  
15cm



**DE-678**  
Hemingway  
15cm



**DE-679**  
Miller  
15cm Str.



**DE-680**  
Miller  
15cm Cvd.



**DE-681**  
Dressing  
Tweezer  
14.5cm



**DE-682**  
College  
Curved  
Tip 15cm



**DE-683**  
Maryam  
Tweezers  
15cm



**DE-684**  
College  
Tweezer  
15cm



**DE-685**  
Perry  
Tweezer  
12.5cm



**DE-686**  
Self-locking  
Tweezer  
15cm

# EXPLORERS


- Зонды
- Sondes
- Zahnsonden
- Зонди
- Sondas
- المستكشفين

Single Ended Probe Solid Handle 4.5mm Length 150mm



Explorers-Single-End


Double-Ended Explorers






- Зонды
  - Sondes
  - Zahnsonden
- Зонди
  - Sondas
  - المستكشفين

Probes Double Ended 6mm, length 170mm



Pointed Probes Double Ended Handle 6mm, length 170mm

**DE-710**  
Williams




# EXPLORERS

- Зонды
  - Sondes
  - Zahnsonden
- Зонди
  - Sondas
  - المستكشفين

EXPLORERS INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

## PLUGGERS




# PERIODONTAL PROBES


- Зонды пародонтологические
  - Sondes parodontales
  - sonda parodontali
- зонди пародонтологічні
  - Sondas parodontales
  - المقاسير اللثة

## SPREADERS

## HEAT CARRIERS



## HEAT CARRIERS




**DE-734**  
ISO  
Size 30  
24mm



- Зонды пародонтологические
- Sondes paradontales
- sonda parodontali
- зонди пародонтологічні
- Sondas paradontales
- المسابير اللثة

## ENDODONTICS INSTRUMENTS

### ANTERIOR



### LUKS





**DE-762**  
Gracey Curettes 1/2



**DE-763**  
Gracey Curettes 3/4



**DE-764**  
Gracey Curettes 5/6



**DE-765**  
Gracey Curettes 7/8



**DE-766**  
Gracey Curettes 9/10



**DE-767**  
Gracey Curettes 11/12



**DE-768**  
Gracey Curettes 13/14



**DE-769**  
Gracey Curettes 1/2



**DE-770**  
Gracey Curettes 3/4



**DE-771**  
Gracey Curettes 5/6



**DE-772**  
Gracey Curettes 7/8



**DE-773**  
Gracey Curettes 9/10



**DE-774**  
Gracey Curettes 11/12



**DE-775**  
Gracey Curettes 13/14



**GRACEY MINI FIVE**

Hollow Handle 9.5mm, Length 170mm Ceramic Coated Handle



**DE-777**  
Fig. 3/4 Gracey 3M



**DE-778**  
Fig. 5/6 Gracey 3M



**DE-779**  
Fig. 7/8 Gracey 3M  
Buccal and lingual/  
Palatal surfaces



**DE-780**  
Fig. 11/12 Gracey 3M  
Mesial



**DE-781**  
Fig. 13/14 Gracey 3M  
Distal

**DE-776**  
Fig. 1/2  
Gracey 3m



**GRACEY AFTER FIVE**

Hollow Handle 9.5mm, Length 170mm Ceramic Coated Handle



**DE-783**  
Gracey 3/4 AF



**DE-784**  
Gracey 5/6 AF



**DE-785**  
Gracey 7/8 AF



**DE-786**  
Gracey 11/12 AF



**DE-787**  
Gracey 13/14 AF

**DE-782**  
Gracey  
1/2 AF





**Gracey Currettes**



$\frac{3...1}{3...1} \mid \frac{1...3}{1...3}$  **DE-788 Fig. 1/2 GR - Gracey**  
Buccal surfaces



$\frac{5...1}{5...1} \mid \frac{1...5}{1...5}$  **DE-789 Fig. 3/4 GR - Gracey**  
Lingual/Palatal surfaces



$\frac{5...1}{5...1} \mid \frac{1...5}{1...5}$  **DE-790 Fig. 5/6 GR - Gracey**



$\frac{8...4}{8...4} \mid \frac{4...8}{4...8}$  **DE-791 Fig. 7/8 GR - Gracey**  
Buccal and lingual/  
palatal surfaces



$\frac{8...4}{8...4} \mid \frac{4...8}{4...8}$  **DE-792 Fig. 9/10 GR - Gracey**  
Buccal and lingual/  
palatal surfaces more angled



$\frac{8...4}{8...4} \mid \frac{4...8}{4...8}$  **DE-793 Fig. 11/12 GR - Gracey**  
Mesial






$\frac{8...4}{8...4} \mid \frac{4...8}{4...8}$  **DE-794 Fig. 13/14 GR - Gracey**  
Distal



$\frac{8...4}{8...4} \mid \frac{4...8}{4...8}$  **DE-795 Fig. 15/16 GR - Gracey**  
Mesial



$\frac{8...4}{8...4} \mid \frac{4...8}{4...8}$  **DE-796 Fig. 17/18 GR - Gracey**  
Distal



**Sickle Scalers**



DE-797 Fig. U15/33 - Towner-Jaquette



DE-798 Fig. U15/30 - Towner-Jaquette



DE-799 Fig. H5/33



DE-800 Fig. H6/H7



DE-801 Fig. 30/33 - Jaquette



DE-802 Fig. T2/T3 - Taylor



DE-803 Fig. 204



DE-804 Fig. 204S




DE-805 Fig. 11/12 - MC - McCall  
Mesial




DE-806 Fig. 13/14 - MC - McCall  
Distal



**Straight and curved sickle scalers**  
Feature: 2 cutting edges curved facial component tipped ends



**Universal curettes**  
are characterized by an angle of 90° between the terminal stem and the cutting component.



**Gracey curettes**  
are characterized by an angle of 70° between the terminal stem and the cutting component.

Universal Curettes



DE-807 Fig. 13S/14S - MC - McCall



DE-808 Fig. 17S/18S - MC - McCall  
Distal



DE-809 Fig. M23



DE-810 Fig. 204SD



DE-811 Fig. YG7/8-Younger Good



DE-812 Fig. L 1/2 - Langer



DE-813 Fig. L 3/4 - Langer



DE-814 Fig. L 5/6 - Langer



DE-815 Fig. M 23A




DE-816 Fig. M 23B




**Straight and curved sickle scalers**

Feature: 2 cutting edges curved facial component tipped ends



**Universal curettes**

are characterized by an angle of 90° between the terminal stem and the cutting component.



**Gracey curettes**

are characterized by an angle of 70° between the terminal stem and the cutting component.



Periodontal knives





Chisels  
Diamond coated margin trimmers

Thanks to the diamond-coated hand instruments the margins are evenly levelled.

Margin trimmers



**DE-835**  
1TH2 Hirt



**DE-836**  
5TH6 Hirt



**DE-837**  
3TH4 Hirt



**DE-838**  
7TH8 Hirt



**DE-839**  
62NS61

2.0 MM



**DE-840**  
37NS36

1.5 MM



**DE-841**  
66NS65

1.0 MM



**DE-842**  
79CI180 Black



**DE-843**  
77CI78 Black



**DE-844**  
279CI280 Black



**DE-845**  
277CI278 Black



**DE-846**  
179CI180 Black



**DE-847**  
177CI178 Black



**CHISELS**



**DE-848**  
1/2 Werdelstaedt



**DE-849**  
3/4 Werdelstaedt



**DE-850**  
7 Straight



**DE-851**  
5/6 Straight



**DE-852**  
7/10 Straight



**DE-853**  
8/9 Binangle



**DE-854**  
11/12 Binangle



**DE-855**  
40/41 Binangle

**Tartar Remover Instruments Solid Handle Ø 6 mm**



**DE-856**  
Fig. 1



**DE-857**  
Fig. 1L  
Left



**DE-858**  
Fig. 1R  
Right



**DE-859**  
Fig. 2L  
left



**DE-860**  
Fig. 2R  
Right



**DE-861**



**DE-862**



**DE-863**



**DE-864**  
WHITE 28



**DE-865**  
WHITE 29



**DE-866**  
WHITE 31



**DE-867**  
WHITE 32



**DE-868**  
WHITE 33



**DE-869**  
WHITE 34



**DE-870**  
WESTCOTT 21



**DE-871**  
WESTCOTT 21B



**DE-872**



**DE-873**



**DE-874**  
WHITE



**DE-875**  
WHITE



DE-876



DE-877



DE-878



DE-879



DE-880



DE-881



DE-882



DE-883



DE-884



DE-885  
PLANATOR



DE-886  
PLANATOR



DE-887



DE-888



DE-889



DE-890



DE-891



DE-892



DE-893



DE-894






**Ball Burnishers**



**Hatches**



## Spatulas rounded and sharp




## Spatulas Blade Thickness 0.22 mm





Spatulas



Spatulas-Pluggers





Water drop Our excavators have particularly sharp blades, free of burrs and perfectly rounded.
























**DE-1005**  
WILLIAMS



**DE-1006**  
WILLIAMS



**DE-1007**  
WILLIAMS



**DE-1008**



**DE-1009**



**DE-1010**



**DE-1011**



**DE-1012**



**DE-1013**



**DE-1014**









**Serrated Pluggers**



|          |                         |          |
|----------|-------------------------|----------|
| Ø 0.9 mm |                         | Ø 1.2 mm |
|          | <b>DE-1034</b><br>OP36M |          |
| Ø 1.1 mm |                         | Ø 1.6 mm |
|          | <b>DE-1035</b><br>OP36D |          |
| Ø 1.6 mm |                         | Ø 2.2 mm |
|          | <b>DE-1036</b><br>Op36  |          |
| Ø 1.0 mm |                         | Ø 1.7 mm |
|          | <b>DE-1037</b><br>Op52  |          |
| Ø 2.2 mm |                         | Ø 2.6 mm |
|          | <b>DE-1038</b><br>Zi51  |          |
| Ø 1.8 mm |                         | Ø 2.2 mm |
|          | <b>DE-1039</b><br>Op44  |          |
| Ø 1.7 mm |                         | Ø 2.6 mm |
|          | <b>DE-1040</b><br>1BR2  |          |
| Ø 2.6 mm |                         | Ø 3.6 mm |
|          | <b>DE-1041</b><br>2BR3  |          |





DE-1042 MF2C  
Fuzzi



DE-1043  
Pj4



DE-1044 MB3C  
Besek



DE-1045  
HZ3C



DE-1046  
CPRO3C



DE-1047 Mo1  
Monari



DE-1048 TH28C  
Hirt



DE-1049  
HZ25C



DE-1050  
Ck10



DE-1051  
CK28BC



DE-1052 Op37  
Carver



DE-1053  
28BR29



### Round excavators

**Cariology**

- for the removal of softened, carious dentin
- for the removal of temporary fillings

**Surgery**

- to sever the periodontal ligament prior to extraction
- for the removal of granulation tissue after extraction

Design: Round working end Curved cutting edge



Round Excavator, 1.0 mm  
DE-1054



Round Excavator, 1.5 mm  
DE-1055



Round Excavator, 2.0 mm  
DE-1056

### Triple-angled back action excavators



Triple-angled Back Action Excavator, 1.0 mm  
DE-1057



Triple-angled Back Action Excavator, 1.5 mm  
DE-1058



Triple-angled Back Action Excavator, 2.0 mm  
DE-1059



Triple-angled Back Action Excavator, 3.5 mm  
DE-1060





# FILLING INSTRUMENTS


- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

COMPOSITE/PLASTIC FILLING INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM

BLACK


PLUGGERS/CONDENSERS

Hollenback



MARQUETTE CONDENSER



MORTONSON



- Инструменты для реставрации
- Instruments pour plomage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

## RESTORATIVE INSTRUMENTS SOLID HANDLE 6MM, LENGTH 170MM Back Action Filling Instruments

### Tanner Filling Instruments



### BALL



### BURNISHERS

### BEAVERTAIL



### ACORN




### LADMORE



# CARVERS INSTRUMENTS

- Инструменты для реставрации
- Instruments a modeler
- Modellierinstrumente
- Инструменти для реставрації
- Instrumentos para modelar
- أدوات تقطيع

## P.K. THOMS WAXING INSTRUMENTS SOLID HANDLE 6mm




DE-1095  
Composite Instruments Blue  
Titanium with Cassettes



DE-1094  
P.k. Thomas Waxing Instruments Kit

## COMPOSITE INSTRUMENTS TITANIUM BLUE ENDODONTIC



## CARVERS INSTRUMENTS

- Инструменты для реставрации
- Instruments a modeler
- Modellierinstrumente
- Инструменти для реставрації
- Instrumentos para modelar
- أدوات تطبيع



DE-1102



DE-1103



DE-1104



DE-1105



DE-1106



Cassette For 5 Instruments  
DE-1107




# FILLING INSTRUMENTS

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

## FILLING INSTRUMENTS (Gold Pointed) SOLID HANDLE 6MM, LENGTH 170MM



- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)
- Инструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك



DE-1118

DE-1119

DE-1120

DE-1121

DE-1122

DE-1123

DE-1124

DE-1125

DE-1126

- Инструменты для реставрации
- Instruments pour plombage
- Füllungsinstrumente (mit Kontrawinkel)

- Інструменти для реставрації
- Instrumentos para obturaciones
- ملء الصكوك

## Composite Gold Tip Instruments





- Прочие товары
- Divers produits et services
- Sonstige Produkte
- Інші товари
- Varios Productos
- المنتجات المتنوعة



**DE-1138**  
X Ray Film Hangers



**DE-1139**  
X Ray Film Hangers



**DE-1140**  
X Ray Film Hangers



**DE-1141**  
X Ray Film Hangers



**DE-1142**  
X Ray Film Hangers



**DE-1143**  
X Ray Film Hangers



**DE-1144**  
X Ray Film Hangers



**DE-1145**



**DE-1146**



**DE-1147**



Dissecting Set  
**DE-1148** Set of 9 Pcs.  
**DE-1149** Set of 11 Pcs.



**DE-1150**  
Napkin Holder with Plastic Chain

**DE-1151**  
Napkin holder with metal chain



Mouth Props  
**DE-1152**  
Small

Mouth Props  
**DE-1153**  
Medium

Mouth Props  
**DE-1154**  
Large

### Key Rings



**DE-1155**



**DE-1156**



**DE-1157**



**DE-1158**




**DE-1159**



**DE-1160**

# CEMENT SPATULAS

- Шпатели для цемента
- Spatules à ciment
- Zementspatel
- Шпатели для цементу
- Espátulas para cemento
- ملاعق الاسمنت



- Ножи для воска
  - Couteaux a cire
  - Wachsmesser
- Ножи для воску
  - Cuchillos para cera
  - السكاكين الشمع

ALGINATE, PLASTER SPATULAS, WAX KNIVES PLASTER FOR KNIVES & PLASTER SHEAR



**DE-1183**  
Wooden Handle  
Fig. 5  
21.5 cm



**DE-1184**  
Wooden Handle  
Fig. 7  
19.5 cm curved



**DE-1185**  
Wooden Handle  
Fig. 8  
20 cm curved



**DE-1186**  
Wooden Handle  
Fig. 9  
22.5 cm curved



**DE-1187**  
Hook Knives



**DE-1188**  
Fig. 1  
Lessmann  
12.5 cm



**DE-1189**  
Fig. 2  
Lessmann  
17 cm



**DE-1190**  
Fig. 1  
with small spoon  
12.5 cm



**DE-1191**  
Fig. 2  
with large spoon  
17 cm



**DE-1192**  
Bacolite Handle  
Wax Knives,  
Small



**DE-1193**  
Bacolite Handle  
Wax Knives,  
Large



**DE-1194**  
Wooden Handle  
Wax Knives,  
Small



**DE-1195**  
Wooden Handle  
Wax Knives,  
Large



**DE-1196**  
Plaster Knife



**DE-1197**  
Fig. 1  
17 cm



**DE-1198**  
Fig. 3  
16 cm



**DE-1199**  
Fig. 4  
16 cm



**DE-1200**  
Special Inox  
19.5 cm



# ARTICULATORS

- Артикуляторы
- Artocilar
- Artikulations

- Артикулятори
- Articuler
- آرتیکولاتورس



**DE-1201**

**Technical Datas:**

|                                  |              |
|----------------------------------|--------------|
| Intercondylar distance           | 120 mm.      |
| condylar inclination             | fixed at 30° |
| Bennet angle                     | fixed at 7°  |
| Incisal table inclination        | 10°          |
| Distance between mounting plates | 85 mm.       |
| Weight                           | 400 gr.      |
| Opening (item 5000)              | 125°         |
| Opening (item 5010)              | 180°         |



**DE-1202**

**Technical Datas:**

|                                  |              |
|----------------------------------|--------------|
| Intercondylar distance           | 110 mm.      |
| condylar inclination             | fixed at 30° |
| Bennet angle                     | fixed at 7°  |
| Incisal table inclination        | flat         |
| Distance between mounting plates | 85 mm.       |
| Weight                           | 450 gr.      |
| Opening                          | 180°         |



**DE-1203**

Stainless Steel Finishing

**DE-1204**

Royal Blue



**DE-1205**



**DE-1206**



**DE-1207**



**DE-1208**



**DE-1209**



**DE-1210**



**DE-1211**




**DE-1212**



**DE-1213**



**DE-1214**  
Nakawalni  
Dental Laboratory  
Instrument



**DE-1215** Adult Palatal /  
**DE-1216** Child Palatal

**DE-1217**  
Buccal

**DE-1218** Adult Palatal /  
**DE-1219** Child Palatal

**DE-1220** Adult Palatal /  
**DE-1221** Adult Buccal




**DE-1222** Lingual upper /  
**DE-1223** Lingual Lower

**DE-1224** Adult Buccal /  
**DE-1225** Child Buccal

**DE-1226** Child Palatal /  
**DE-1227** Child Buccal

**DE-1228**  
Palatal-Adult



**DE-1229**  
Palatal-Adult Large

**DE-1230**  
Palatal-Child

**DE-1231**  
Buccal-Small

**DE-1232**  
Buccal

**DE-1233**  
Lingual

**DE-1234**  
Lingual-Large